SAPPHIRE

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SAPPHIRE

Light Cure Nano - Composite

Sapphire is a **visible light cure nano - composite** for **anterior and posterior** restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-0.9 µm. Filler content: 81 % (wt), 63 % (vol.). The composit material meets the requirements of **DIN EN ISO 4049**.

Indications:

Suitable for all cavities.

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of hard setting calciumhydroxid material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.

Preparation of Tooth and Colour Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and colour matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate **Sapphire** shade.

Sapphire is shaded according to VITA®-Shades. It is also available in special shades and bleach shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

Best pulp protection in case of deep excavations is the photo-activated disinfection (**Lazurit®**). If this is not possible, deep excavations should be covered with a thin layer of hard setting calciumhydroxid material.

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with **Sapphire Conditioner** and subsequent application of a light cure bonding agent (e.g. **Sapphire Bond**) is recommended. Otherwise use a self etching bonding agent (e.g. **Sapphire Bond SE**) before application of the **Sapphire**.

Application and Curing of SAPHIRE

Apply **Sapphire** in the selected shade. **Sapphire** is easy to model. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite (e.g. **Sapphire Flow**), followed by the application of **Sapphire**.

In case of small cavities **Sapphire** can be applied directly.

For an optimum result apply **Sapphire** in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity of **at least 1000mW/cm²**, cure each increment as follows:**

Lighter shades (e.g. А1, А2, С2) **20 seconds**

Darker shades (e.g. А3.5, В3) **30 seconds**

Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

Storage

Do not store above 25 °C (78 °F). Cool storage prolongs shelf life.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.

- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Distributor:	TBI® Company, 6/2 Ostozhenka Str. RU-119034 Moscow	
Manufacturer:	S&C Polymer GmbH, Elmshorn, Germany	